

ACTREC

TEST DIRECTORY REPORT

LOCATION DETAILS

Department name : CANCER CYTOGENETICS DEPARTMENT
Location : Room No. 6, Ground Floor, CCE Building.

FUNCTIONAL DETAILS

Test Name : XX/XY (Chimerism Studies) in Sex mismatch Bone Marrow Transplantation (BMT)
Test code : T470
Test status : ACTIVE
Modality : FI
Quantity required : Bone marrow: 1-4 ml, Peripheral Blood: 4-5 ml
Collection Instructions : Refer Primary sample collection manual
Nature of specimen : Bone marrow, Peripheral Blood
Method used : Fluorescence in situ hybridization (FISH)
Type of container : Sterile sodium heparin Green top vacutainer
Transport instructions : Refer Primary sample collection manual
Time for additional exam on stored sample, if applicable : NA
Resource person for report status : Officer-in-charge
Resource person for test query : Officer-in-charge
Patient instruction : Refer Primary sample collection manual
Dept. acceptance time : Monday to Friday-9:30 a.m. till 5:00 p.m.
Saturday- 9:30 a.m. till 12:00 p.m.
Turn around time : 1-2 working days
Sample storage after reporting : Fixed pellet stored in minus 20 deep freezer
Telephone : 022-27405000
Test schedule : Daily
Extension : 5758, 5759
Use : Sex mismatch chimerism